ABSTRACT

# CHARACTERIZATION OF TP53 MUTATIONS BY DNA AND RNA SEQUENCING OF HIGH GRADE SEROUS OVARIAN CANCER (HGSOC) FFPE SECTIONS

<u>Diane Ilsley</u><sup>1</sup>, Marie Fahey<sup>1</sup>, Brian C Haynes<sup>1</sup>, Gary J Latham<sup>1</sup>, Huiping Zhu<sup>1</sup>, Andrew Kincaid<sup>1</sup>, Jeffrey Shelton<sup>1</sup>, Elizabeth B Somers<sup>2</sup>, and Daniel J O'Shannessy<sup>2</sup> Asuragen, Inc., Austin, TX; <sup>2</sup>Morphotek, Inc., Exton, PA

## SUMMARY


- Primary FFPE tumor samples from a cohort of >200 high grade serous ovarian cancer (HGSOC) patients were profiled using an integrated analysis comprised of both DNA and RNA sequencing to gain a more complete understanding of the functional consequences of *TP53* mutations.
- TP53 mutation status was determined using the QuantideX® NGS TP53 Coding Exon Panel and the AmpliSeq® Cancer Hotspot Panel, and regions of TP53 covered by both assays were found to be concordant.
- Mutation status determined using matched RNA-Seq data was found to be highly concordant with the QuantideX® NGS panel.
- RNA-Seq gene expression data was integrated with the DNA mutation data resulting in the identification of a set of high confidence associations between *TP53* mutations and RNA expression changes in genes regulated by p53.

### INTRODUCTION

TP53 is the most frequently mutated gene across all cancers at approximately 50%, and in cancers such as HGSOC, up to 96%. In addition to the high incidence of DNA mutations at the gene locus, the TP53 pathway is also dysregulated through alternative mechanisms which involve other members of the regulatory cascade, such as overexpression of MDM2. Asuragen has developed sensitive NGS assays and algorithms to characterize DNA mutations and RNA expression profiles from the most challenging clinical samples, such as archival FFPE samples. In order to enable a more complete understanding of the functional consequences of TP53 mutations (or lack thereof), we present an application of these approaches to characterize TP53 function in FFPE sections from HGSOC through an integrated analysis of DNA and RNA sequencing.

#### METHODS

Primary debulking surgery FFPE tumor sections were collected from a cohort of 235 HGSOC patients. DNA from macro-dissected tumor and matched normal blood specimens were profiled by the QuantideX® NGS *TP53* Coding Exon Assay (Asuragen, Inc.), and library analysis was performed using the QuantideX® NGS Reporter. *TP53* mutation status of the tumor specimens was independently assessed by the AmpliSeq Cancer Hotspot Panel v1 (Thermo Fisher Scientific, Inc.). Whole transcriptome RNA-Seq was performed using ribosomal RNA-depleted total RNA, strand-specific libraries, and 50 bp paired-end sequencing. Gene and isoform expression quantitation as well as SNV and indel calling was performed on the RNA-Seq libraries.



**Figure 1. QuantideX and AmpliSeq NGS panel gene coverage of TP53.** The coding exons are represented in blue and the noncoding exons in orange. The QuantideX TP53 panel covers 100% of the TP53 coding bases, 99.5% of all COSMIC variants and all ovarian TCGA TP53 variants (390). In contrast, the AmpliSeq Hotspot panel covers 68.4% of the coding bases, 68.4% of the TP53 COSMIC variants and 303 ovarian TCGA variants.

#### RESULTS




Figure 2. Variants identified by the QuantideX NGS panel aligned to the p53 protein domains. The plot shows the distribution of both SNVs (187 somatic) and indels (37), which were found to be mutually exclusive at the patient level. Consistent with previous studies, TP53 mutations were present in 84% of the HGSOC cohort. Symbol color indicates the type of mutation and symbol shape indicates variant annotation.





Figure 3. Positive call summary comparing results by panel. A) The Venn diagram shows the concordance of variant calls in gene regions covered by both QuantideX and AmpliSeq after filtering for germline variants and applying a coverage filter (100X). RNA-Seq data confirmed 6/26 positive calls in QuantideX and 2/31 positive calls in AmpliSeq. The variants called positive by AmpliSeq (29/31) were not called positive by QuantideX because they did not meet one of the following criteria: sensitivity threshold, coverage filter, or contain matching forward and reverse reads. Considering the total content coverage of the QuantideX Panel, QuantideX called an additional 51 variants that were not covered by the AmpliSeq panel. B) Variant frequency comparison of QuantideX (Y-axis) versus AmpliSeq (X-axis) for the 110 variants called by both panels. Variants are colored by mutation functional class.



**Figure 4. QuantideX identifies positive variants missed by AmpliSeq.** The variant frequency by QuantideX (Y-axis) is plotted versus the variant frequency by AmpliSeq (X-axis) highlighting the variant calls missed by AmpliSeq due low sequencing coverage (very low = <100X) or below the sensitivity threshold. The symbol shape is variant annotation and symbol color reflects depth of sequencing coverage in AmpliSeq. The top panels show the variants that were confirmed by RNA-Seq while the bottom panel contains variants not confirmed by RNA-Seq due to low sequencing coverage (2-20X).




Figure 5. TP53 variant type differentially affects TP53 expression and expression of TP53 regulated genes. The graph shows boxplots of the median log2 CPM for a subset of genes (indicated along the X-axis) directly or indirectly regulated by p53. TP53 variant type is indicated by color and boxplots shaded highlight the genes and TP53 variants that are differentially expressed relative to wild-type TP53 based on a FDR corrected p-value < 0.05.

## CONCLUSIONS

- We successfully used NGS to characterize the DNA mutations and RNA expression profiles from a cohort of 235 archival HGSOC FFPE specimens.
- Because *TP53* is highly mutated across all exons, a hot spot mutation panel does not provide adequate gene coverage for mutation assessment. The QuantideX® NGS *TP53* Panel identified mutations in 60 additional patients that would have been missed by a hotspot panel.
- Integration of *TP53* DNA mutation data with matching RNA-Seq data provides a more complete understanding of the functional consequences of the different classes of *TP53* mutations through association of gene signatures with gain and loss of p53 function.



